**ORDERING INFORMATION:**

To place an order for URINE CREATININE CALIBRATOR, contact the Customer Service Department at 1 (800) 526-3821.

Catalog No.: DR0091  
 Description: URINE CREATININE CALIBRATOR  
 Packaging: 120 mL

**ASSAY VALUES**

LOT 2109038A EXP. DATE: 10/31/23

Beckman Coulter AU® Series Analyzers

Test Name	Constituents	Units
CRE1U	Creatinine	umol/L 8354 mg/dL 94.5

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
**GHS Hazard Classification**

DANGER

H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves, protective clothing and eye/face protection.
P301+P330+P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303+P361+P353	IF ON SKIN (or hair): Rinse skin with water.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/physician.
Hydrochloric Acid 1-3%	
<b>SDS</b> Safety Data Sheet is available at <a href="http://techdocs.beckmancoulter.com">techdocs.beckmancoulter.com</a>	

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Manufactured for:  
 Beckman Coulter, Inc.  
 250 S. Kraemer Blvd.  
 Brea, CA 92821  
 Made in USA

## INTENDED USE

URINE CREATININE CALIBRATOR is intended for use as a urine chemistry calibrator for methods run on the BECKMAN COULTER AU® series of Clinical Chemistry Systems (refer to "Use Limitations" section below).

## CONSTITUENTS

URINE CREATININE CALIBRATOR is an aqueous based product assayed for Creatinine.

## PRECAUTIONS

- 1 FOR IN VITRO DIAGNOSTIC USE
- 2 DO NOT INGEST-Toxicity has not been established.

## PREPARATION

URINE CREATININE CALIBRATOR is ready to use. No preparation is needed.

## STORAGE AND STABILITY

Store unopened bottles of URINE CREATININE CALIBRATOR at 2 to 8°C. When stored as directed, the calibrator is stable until the expiration date printed on the label.

Once opened, the URINE CREATININE CALIBRATOR is stable for 8 months.

## USE LIMITATIONS

This calibrator has not been tested for use with any other Chemistry Systems other than those listed on the Assay Value section.

For best results avoid prolonged exposure of the samples to air. Withdraw sample and replace cap immediately. Maintain the bottle at 2 to 8°C. Assay the sample without delay.

The results obtained using this calibrator are dependent upon several factors, including proper storage of the calibrator, proper technique, and use of Beckman Coulter AU Clinical Chemistry Analyzers and their respective reagents.

The DxC 700 AU is not available in all geographies.

## INDICATIONS OF DETERIORATION

Bacterial contamination produces an increase in turbidity, a characteristic odor, and/or a decrease in the creatinine level. Discard the bottle if any of these phenomena are observed.

## RECOMMENDED PROCEDURE

- 1 Gently invert the bottle several times to ensure a homogeneous mixture.
- 2 Transfer sufficient volume of the calibrator to sample cups. Handle this calibrator with the same care used for patient samples.
- 3 Replace the cap immediately and store unused calibrator at 2 to 8°C.
- 4 Refer to the appropriate Instrument User's Guide for System Calibration Information.
- 5 Good Quality Control Practices should be observed to assure proper System performance.
- 6 Do not pour previously dispensed calibrator back into original bottle.
- 7 The 2 dimensional bar code is intended for use with the DxC 700 AU analyzer only.

## VALUE ASSIGNMENT

The assay value has been determined in accordance with BECKMAN COULTER test protocol (adopted from CLSI) and it applies only to this lot of calibrator. The creatinine value for each AU Clinical Chemistry Analyzer has been assigned using its respective reagent. Any instrument or reagent modifications may invalidate these assay values.

**PLEASE ENSURE THAT THE LOT NUMBER ON THE CALIBRATOR BOTTLE IS THE SAME AS THE ONE LISTED IN THE ASSAY VALUE SECTION OF THIS PACKAGE INSERT.**

The assay value for the constituent is traceable to the National Institute of Standards and Technology (NIST SRM 3667).

## TECHNICAL ASSISTANCE

For technical assistance contact the BECKMAN COULTER Technical Services Department at 1 (800) 223-0130.

## REFERENCES

1. National Institute of Standards and Technology, Office of Standard Reference Material, Room 8322, Chemistry Building, Gaithersburg, MD 20899.

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